Using BiSEp to nominate candidate Synthetic Lethal gene pairs

Mark Wappett

January 25, 2017

1 Introduction

BiSEp (Bimodality subsetting expression) is a set of computational tools that enable the user to nominate candidate synthetic lethal (SL) gene pairs. The toolkit is based on the assumption that a clear on / off gene expression profile is indicitive of tumour loss, and is detectable as bimodality or non-normality.

This vignette demonstrates how the toolkit can be used to nominate, assess and visualise candidate SL pairs nominated from gene expression and mutation datasets.

2 Importing gene expression data

Processed gene expression data from most platforms can be input. We recommend that values are all distributed above zero and are in the log2 scale. Example input data format is shown below:

All input data should be read in a gene by sample format. Our dataset is gene expression data from the Cancer Cell Line Encyclopedia (CCLE) [1], normalised using fRMA [2] and scaled.

3 Identifying bimodal genes in expression data

We next take the processed data matrix and run the bimodal detection tool across it. This generates a list object containing three matrices, the third of which is called DATA and is simply a capture of the input data matrix.

```
> BISEP_data <- BISEP(INPUT_data)
> biIndex <- BISEP_data$BI
> bisepIndex <- BISEP_data$BISEP</pre>
```

The output data frames called biIndex and bisepIndex are the output from the bimodal index function [3] and the novel BISEP function. The output is displayed below:

> biIndex[1:10,]

```
ΒI
             mu1
                      mu2
                              sigma
                                       delta
                                                      pi
MICB
       0.9321880 3.870030 0.7636438 3.847136 0.28249286 1.7320261
YAP1
       0.3058179 3.189638 0.6001395 4.805250 0.11013333 1.5043116
BOK
       1.3351088 3.742910 0.6752271 3.565913 0.18944756 1.3973518
       0.5245744 3.470745 0.6534969 4.508315 0.18604307 1.7543722
MTAP
EPHB2
       2.1856093 4.163538 0.7103121 2.784591 0.34551628 1.3241743
       3.2415304 4.318537 0.7020134 1.534169 0.76649523 0.6490464
BRCA2
TUSC3
       0.5724434 3.854679 0.6501654 5.048309 0.25574133 2.2024615
PHLDA3 1.9027926 4.050551 0.8063510 2.663553 0.37025203 1.2861553
MLH1
       0.5957125 5.146278 0.4644327 9.798116 0.03846154 1.8842532
BRCA1
       3.2805696 4.988761 0.5504041 3.103522 0.06713879 0.7766944
```

> bisepIndex[1:10,]


	V1	V2
MICB	1.448461	0.49751244
YAP1	1.052045	0.34129693
BOK	1.982416	0.63694268
MTAP	0.966708	0.15243902
EPHB2	2.173364	0.71428571
BRCA2	5.170457	0.94339623
TUSC3	2.363063	0.23752969
PHLDA3	5.379839	1.0000000
MLH1	2.543707	0.07518797
BRCA1	3.609235	0.78125000

The biIndex matrix contains all the bimodal scoring information provided to us by the bimodal index function. This includes the delta (distance between two expression modes), pi (proportion of samples in each expression mode) and BI. When combined, these give us an optimal assessment of bimodality in expression data. The bisepIndex function provides a p-value score for nonnormality (column 2) and accurately pin-points the mid-point between the two expression modes for a gene (column 1).


TUSC3 scores the highest in the biIndex table - the density distribution below highlights this:

> plot(density(INPUT_data["TUSC3",]), main="TUSC3 Density Distribution")

TUSC3 Density Distribution



MLH1 Density Distribution



N = 442 Bandwidth = 0.1372

By comparison, MLH1 does not score high for bimodality - but has the lowest p value for non-normality. The density plot demonstrates the unbalanced nature of, and distance between the two populations in a typical non-normal distribution.

4 BIGEE: Bimodality in Gene Expression Exclusivity

Here we take the bimodal / non-normal output from the BISEP tool, and use it as input to the first of the two candidate synthetic lethal detection tools. There are four sample input options to this tool based on the sample type and sample numbers **cell line, cell line low, patient and patient low**. When sample numbers are below 200 we recommend using the input parameters with the low suffix in order to prevent a high false positive rate.

> BIGEE_out <- BIGEE(BISEP_data, sampleType="cell_line")</pre>

- [1] "Selected CELL LINE sample type"
- [1] "Subsetting bimodal index"
- [1] "Filtering"
- [1] "Setting up synthetic lethal detection"
- [1] "Running SL detection"
- [1] "Number of bimodal genes: 12"

	Ι	0%
 =====	Ι	8%
 ============	Ι	17%
 ====================================	Ι	25%
 ===================================	Ι	33%
 ====================================	I	42%
 ==================================		50%
' ==================================		58%
 	1	
	і	67%
=====================================	I	75%
====================================	Ι	83%

l		l	92%
	 ===================================	I	100%
[1]	"Summarising"		


The percent completion graphic displays the progress of the SL detection component of the tool. This will typically take longer the larger the dataset is, and the more bimodal genes that there are. The output from this tool is a matrix containing gene pairs that look potentially synthetic lethal in the dataset, along with a score.

```
> BIGEE_out[1:4,]
```

	gene2	gene	score
2	MTAP	MLH1	11.686754
1	MLH1	BOK	9.995059
9	YAP1	MLH1	7.353136
7	SMARCA4	SMARCA1	6.814576

It is possible to visualise any candidate relationships using the expressionPlot function:


> expressionPlot(BISEP_data, gene1="SMARCA4", gene2="SMARCA1")



SMARCA4 vs. SMARCA1 Log2 Gene Expression plot

and look for those gene pairings that ideally are never expressed at low levels together - the signature that we propose could be indicative of synthetic lethality.

> expressionPlot(BISEP_data, gene1="MTAP", gene2="MLH1")



MTAP vs. MLH1 Log2 Gene Expression plot

5 BEEM: Bimodal Expression Exclusive with Mutation

Here we take the bimodal / non-normal output from the BISEP tool, and use it as input to a tool that detects mutual exclusive loss between bimodally expressed genes and mutated genes. Again, there are four sample input options to this tool based on the sample type and sample numbers **cell line**, **cell line low**, **patient and patient low**. Additionally we also require a second input matrix containing discreet mutation call information. This matrix must be in the rownames = genes, colnames = samples format and there must be overlap between sample names in this mutation matrix, and sample names in the INPUT data matrix seen earlier. The calls in this matrix must be either WT or MUT as shown below:

```
> data(MUT_data)
> MUT_data[1:4,1:10]
```

	5637	42MGBA	639V	647V	769P	8305C	8505C	8MGBA	A101D	A2058	
BRCA2	WT	WT	WT	WT	WT	WT	WT	WT	WT	WT	
PBRM1	WT	WT	WT	MUT	WT	WT	WT	MUT	WT	WT	
SCN2A	WT	MUT	MUT	WT	WT	WT	WT	WT	WT	WT	
CACNA1D	WT	WT	WT	MUT	WT	WT	WT	WT	WT	WT	

Now we can run the function by doing the following:

> BEEMout <- BEEM(BISEP_data, mutData=MUT_data, sampleType="cell_line", minMut=40) [1] "Minimum number of mutations considered for each gene is: 40" [1] "Selected CELL LINE sample type" [1] "Number of bimodal expression genes : 12" [1] "Number of mutation genes wih frequency greater than 40 : 4" 0% 8% 1 = | 17% ========= 25% L _____ _____ 33% 42% 1 _____ 50% _____ 58% 1 67% 1 75% _____ 83% _____ 92% ______ Ι

[1] "Summarising..."


As with the BIGEE tool, the percent completion graphic displays the progress of the SL detection component of the tool. The output from the tool is a matrix containing the gene pairs that look potentially synthetic lethal, along with a number of other columns of metadata including size of high and low expression population, numbers of those populations that are mutant.

Gene1 Gene2 L			
	owerExpressionMutation	Count HighExpression	MutationCount
2 MICB PBRM1	-	0	42
9 BOK BRCA2		0	55
20 EPHB2 CACNA1D		0	41
7 YAP1 SCN2A		0	59
4 MICB CACNA1D		3	38
12 BOK CACNA1D		2	39
5 YAP1 BRCA2		2	53
35 BRCA1 SCN2A		0	59
13 MTAP BRCA2		5	50
18 EPHB2 PBRM1		3	39
Fishers P Value	Percentage of lower s	amples mutated	
2 0.0002127793		0	
9 0.0002552670		0	
20 0.0018895294		0	
7 0.0019129823		0	
4 0.0576643837	3	.7037037037037	
12 0.0973332166	3.	17460317460317	
5 0.0975156120	4.	444444444444	
35 0.1490033473		0	
13 0.1701101513		944444444444	
18 0.2438573394		76190476190476	
_	gh samples mutated Siz	e of low expression	
2	11.6343490304709		81
9	14.5118733509235		60
20			63
	10.8179419525066		63
7	14.8614609571788		63 45
7 4	14.8614609571788 10.5263157894737		63 45 81
7 4 12	14.8614609571788 10.5263157894737 10.2902374670185		63 45 81 63
7 4 12 5	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844		63 45 81 63 45
7 4 12 5 35	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226		63 45 81 63 45 18
7 4 12 5 35 13	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135		63 45 81 63 45 18 72
7 4 12 5 35 13 18	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185		63 45 81 63 45 18
7 4 12 5 35 13 18 Size of high ex	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En		63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361	HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379	HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 379	HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 379 397	HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7 4	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 379 397 361	HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7 4 12	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 397 361 379	HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7 4 12 5	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 397 361 379 397	HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7 4 12 5 35	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 397 361 379 397 424	HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72
7 4 12 5 35 13 18 Size of high ex 2 9 20 7 4 12 5	14.8614609571788 10.5263157894737 10.2902374670185 13.3501259445844 13.9150943396226 13.5135135135135 10.2902374670185 pression population En 361 379 397 361 379 397	HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched HighEnriched	63 45 81 63 45 18 72


> BEEMout

Gene pairs where the mutant gene2 is exclusively mutated, or significantly enriched for mutation in the high expression mode of expression gene1 are those that we propose as candidate SL pairs. It is another manifestation of the neverlow-together relationship we were looking for in the expression data above. We can visualise these gene pairs using the waterfall plotting function built into the package

> waterfallPlot(BISEP_data, MUT_data, expressionGene="MICB", + mutationGene="PBRM1")



> waterfallPlot(BISEP_data, MUT_data, expressionGene="BOK", + mutationGene="BRCA2")



The left panel is the density distribution of the bimodal / non-normal expression gene. The right hand panel is a bimodal-mid-point-centered barplot coloured by the mutation status of the mutation gene.

6 FURE: Functional redundancy between synthetic lethal genes

It is assumed that either gene in a synthetic lethal pair is able to functionally compensate for the loss of the other. We developed this tool to enable the user

to prioritise gene pairs that have some sort of biological redundancy and score these according to gene ontology [4,5].

The tool takes as input either the output from the BIGEE or the BEEM tools. The following example is run on the first couple of results from the BIGEE output

```
> fOut <- FURE(BIGEE_out[1,], inputType="BIGEE")
> frPairs <- fOut$funcRedundantPairs
> allPairs <- fOut$allPairs</pre>
```

```
> allPairs[1,]
```

```
gene2genescoreredundant_idsredundant_termsMolecularFunctionScore2MTAPMLH111.68675GD:0005634nucleus0.13BiologicalProcessScoreCellularComponentScore20.3860.365
```

7 References

1. Berretina J et al. (2012) The Cancer Cell Line Encyclopedia enables predictive modelling of anticancer drug sensitivity. Nature, 283:603-607.

2. McCall MN, Bolstad BM, and Irizarry RA. (2010) Frozen Robust Multiarray Analysis. Biostatistics, 11(2):242-253.

3. Wang J et al. (2009) The Bimodality Index: A criterion for discovering and ranking bimodal signatures from cancer gene expression profiling data. Cancer Informatics, 7:199-216.

4. Carlson, M. (2013) Org.Hs.eg.db: Genome wide annotation for human. R package version 2.8.0.

5. Guangchuang, Y et al. (2010) An R packahe for measuring semanic similarity among GO terms and gene products. Bioinformatics, 26(7), 976-978.